



SPECIFICATION





Mass Product

Contents

1. Description	
1.1 General Description	
1.2 Features	
1.3 Application	
1.4 Features Introduction	
2. Product Specification	
2.1 Size chart	
2.1 Product Parameters	
.application	6
Packing Criterion	7
5. Handling Precautions	8
5.1 Product storage and maintenance	8
5.2 Handling Precautions	9








2. Product Specification


2.1 Size chart




size : 280*208*110mm

2.1 Product Parameters

Product Parameters		
1	Refond PN	RF-MTQAB7N03-B1
2	Battery	3.7V/2300mAh



application



Mode 2: Blu light



Mode 3: red light



Packing Criterion

4.1 Package Diagram



Fig.3-2 Title



Fig.3-2 Title

5. Handling Precautions

5.1 Product storage and maintenance

Do not expose the product to rain, humidity and high temperature.

Open flame sources, such as lit candles, should not be placed on the device

Careful supervision is required to prevent children from using this product.

If not used for a long time, charge the battery at least once every 3 weeks to maintain the battery life

Do not poke this product with sharp metal objects such as keys, hairpins, necklaces or nails.

5.2 Handling Precautions

Products for facial care, not recommended for body use.

Do not use this device for skin allergies, wounds, severe acne, burns, and bacterial skin infections.

Please consult a dermatologist before using this device if you are being treated for a skin disease or have a skin condition.

Not recommended for pregnant, menstruating, lactating women (unstable hormonal balance may cause temporary skin problems).

Do not use for people who are unable to express themselves clearly (including babies).

People with light sensitivity should not use this product if they are allergic to ultraviolet light.

Do not use this product in the rain, damp (bathtub, shower, sink), high temperature (open flame source, oven) environment.

Do not disassemble without permission, there is no user repairable parts inside, related maintenance, please consult qualified maintenance personnel.

Version History/

Date	Revisor	Version	Verifier	Remarks
2021-8-11		E/0		The first edition



Declare

This specification is written both in English and in Chinese and the latter is formal.